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
INDEX	1_ INTRODUCTION	22
	2_ SPECIFICATIONS	22
	2.1_ Serie BAR	22
	2.2_ Common specifications	22
	3_ INSTALLATION	23
	3.1_ General installation	23
	3.2_ Installation of BAR series	23
	4_ KEY DESCRIPTIONS	23
	5_ OPERATION	24
	5.1_ Zeroing the display	24
	5.2_ Taring	24
	5.3_ Weighing a sample	24
	5.4_ Percent weighing	24
	5.5_ Parts counting	25
	5.6_ Accumulated	25
	6_ BATTERY OPERATION	25
	7_ INTERFACE	26
	8_ PARAMETERS	27
	9_ APPENDIX	28
	9.1_ Error code	28
	9.2_ Weighing unit listing	29
	10_ TROUBLE SELF CHECKING	29
	11_ WARRANTY	30

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1. INTRODUCTION


The BAR series of electronic balances provides an accurate, fast and versatile series of general purpose balances with counting and % weighing functions.

There are 2 or 4 models in this serie, with capacities up to 6000g.

They all have stainless steel weighing platforms on an ABS base assembly.

All the keypads are light touch switches and the displays are large easy to read liquid crystal type displays (LCD). The LCD's are supplied with a backlight and have loading parts counting.

All units include automatic zero tracking, automatic tare, and an accumulation facility that allows the count to be stored and recalled as an accumulated total.

2. SPECIFICATIONS

2.1 SERIE BAR

N° Model	SERIE BAR						
	BAR300	BAR600	BAR1500	BAR3000	BAR6000	BAR150I	BAR1500I
Capacity (g)	300	600	1500	3000	6000	150	1500
Readability (g)	0,005	0,010	0,020	0,050	0,100	0,001	0,010
Resolution	1/60.000	1/60.000	1/60.000	1/60.000	1/60.000	1/150.000	1/150.000
Pan size	Ø 120 mm		140 x 150 mm			Ø 80 mm	Ø 120 mm
Wind protector	Yes	Yes	-	-	-	Yes	Yes
Units of measurement	g / ct / lb / oz / d / gn / ozt / dwt / mm / tl. T / tl. c / tl. t / t.						

2.2 COMMON SPECIFICATIONS

Display	16'5 mm LCD, with white color LED backlight
Operate temperature	0°C ~ + 40°C
ADC	Σ-Δ
ADC update	≤ 1/10 second

Housing	ABS plastic
Max. Divisions	1/60.000, 1/150.000
Power	AC adapter (9V/800 mA) Rechargeable battery (6V/1.2Ah)
Gross weight	≤ 1500 = 2'9 kg ≥ 1500 = 2'7 kg

3. INSTALLATION

3.1. GENERAL INSTALLATION

The scales should be sited in a location that will not degrade the accuracy.

Avoid extremes of temperature. Do not place in a direct sunlight or near air conditioning vents.

Avoid unsuitable tables. The tables or floor must be rigid and not vibrate. Do not place near vibrating machinery.

Avoid unstable power sources. Do not use near large users of electricity such as welding equipment or large motors.

Avoid high humidity that might cause condensation. Avoid direct contact with water. Do not spray or immerse the scales in water.

Avoid air movement such as from fans or opening doors. Do not place near open windows.

Keep the scales clean.

Do not stack material on the scales when they are not in use.

3.2. INSTALLATION OF BAR SERIES

Level the scale by adjusting the four feet. The scale should be adjusted such that the bubble in the spirit level is in the center of the level and the scale is supported by all four feet. If the scale rocks readjust the feet.

Attach the adapter to the connector on the side of the scale. Please use the power as the label.

For balances with wind shielding, please put wind shielding on the top of balance, just check reading after load weight and close wind shielding.

4. KEYPAD DESCRIPTION

KEY	DESCRIPTION	
	PRIMARY FUNCTION	SECONDARY FUNCTION
	Set the zero point for all subsequent weighing. The display shows zero.	“Enter” key when setting parameters or other functions.
	Tares the scale. Stores the current weight in memory as a tare value, subtracts the tare value from the weight and shows the results. This is the net weight.	

ENG FR ES		Enters the percent weighing function. When in percentage mode and current value is zero, press this key to return normal weighing mode.	
		To print the results to a PC or printer using the optional RS-232 interface. It also adds the value to the accumulation memory if the accumulation function is not automatic.	"ESC" is to return to normal operation when the scale is in a parameter setting mode.
		This key will select either kilograms, pounds, ounce for the weighing unit when weighing mode. When counting mode, allows the weight, unit weight, and count to be seen when parts counting.	A secondary function, change current value for parameters or other functions.
		Use this key to enter counting mode when normal weighing mode.	Use this key to return normal weighing mode when counting mode.
		Use this key to turn on or turn off the power.	

5. OPERATION

5.1. ZEROING THE DISPLAY

You can press the  key at any time to set the zero point from which all other weighing and counting is measured, within 10% of power up zero. This will usually only be necessary when the platform is empty. When the zero point is obtained the display will show the indicator for zero.

The scale has an automatic rezeroing function to account for minor drifting or accumulation of material on the platform. When the weight goes over 20 % of the full capacity, "read zero" function is obsolete.

5.2. TARING

Zero the scale by pressing the  key if necessary. The zero indicator will be on.

Place a container on the platform, a value for its weight will be displayed.

Press the  key to tare the scale. The weight that was displayed is stored as the tare value and that value is subtracted from the display, leaving zero on the display. The "TARE" indicator will be on. As product is added only the weight of the product will be shown. The scale could be tared a second time if another type of product was to be added to the first one. Again only the weight that is added after taring will be displayed.

When the container is removed a negative value will be shown. If the scale was tared just before removing the container this value is the gross weight of the container plus all product that was removed. The zero indicator will also be on because the platform is back to the same condition it was when the  key was last pressed.

5.3. WEIGHING A SAMPLE

To determine the weight of a sample first tare the empty container then place the sample in the container. The display will show the weight and the units of weight currently in use.

5.4. PERCENT WEIGHING

The scale will allow a sample weight to be shown as 100%. Then any other weight placed on the scale will be displayed as a percentage of the original sample. For example is 350g is placed on the scale and the  key is pressed the display will show 100.00%.

Removing the 350g weight and putting a 300g weight on the scale the display will show 85.71% as 300g is 85.71% of 350g.

Note: the scale may jump by large numbers unexpectedly if small weights are used to set the 100% level. For example if only 23.5g is on a scale with 0.5g increments and the scale is set to 100%, the display will show 100.00%, however a small change of weight will cause the display to jump to 102.13% as one scale division (0.5g) increase to 24.0g will be equivalent to a 2.13% increase.

Removing the 350g weight and putting a 300g weight on the scale the display will show 85.71% as 300g is 85.71% of 350g.

Pressing the  key will return the scale to weighing.

5.5. PARTS COUNTING

When the scale is showing weight, pressing the  key will start the parts counting function.

Before beginning, tare the weight of any container that will be used, leaving the empty container on the scale. Place the number of samples on the scale. The number should match the options for parts counting, 10, 20, 50, 100 or 200 pieces.

Press the SMPL key to begin. The scale will show "SP 10" asking for a sample size of 10 parts. You can press  key to select sample quantity: 10, 20, 50, 100, 200 and back to 10.

Press the  when the number matches the number of parts used for the sample. As more weight is added the display will show the number of parts (pcs).

Press the  key to display unit weight (g/pcs), Total weight (kg) or the count (pcs). Press the  key to return to normal weighing.

5.6. ACCUMULATION

When the scale is set to manual accumulation the weight displayed will be stored in memory when the  key is pressed and the weight is stable.

The display will show "ACC 1" and then the total in memory for 2 seconds before returning to normal. If the optional RS232 interface is installed the weight will be output to a printer or PC.

Remove the weight, allowing the scale to return to zero and put a second weight on. Press , the display will show "ACC 2" and then show the new total. Continue until all weights have been added.

To view the totals in memory press the  key when the scale is at zero. The display will show the total number of items "ACC xx" and the total weight before returning to zero. The totals will also be printed via the RS232 interface.

To erase the memory press  to view the totals and then press the  key to clear the memory.

* When you don't connect to printer, please don't set RS232 to LP50 mode, otherwise accumulate operate will be slow. See detail in section 8.

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The balance come with a rechargeable battery, can be operated from the battery if desired. The battery life is approximately 40 hours.

There is a battery symbol on display used to indicator battery status:

Flick: charging battery

Full: battery full or charge finished

Half: battery with half power

Blank: need be charged

When the battery needs to be recharged, the indicator lights. Please, load the battery or use an adapter of alternating current when he/she sees this lit indicator.

The scales can work with batteries if one wants. The duration of the battery is of approximately 72 hours.

To load the battery plug it to the electric net. It is not necessary to light the scale. The recharge should last 12 hours to fill the battery to full capacity.

There is a LED that indicates the state of the load of the battery. When the scale is plugged to the electric net, the battery goes into it will load. If LED is green the battery has complete load. If it is red, the battery needs load, and the yellow indicates that the battery is in recharge process.

With the use, the battery can go losing capacity to make complete loads. If the duration of the load of the battery becomes unacceptable, contact its distributor.

Note: the new batteries are partially loaded. Before being able to use their scale, it should install and to load the battery following the following instructions.

Some batteries have a better yield after several complete cycles of load / it discharges.

The yield of the batteries depends on numerous factors, included the configuration of the backlight and the operative one.

Not never use a loader or a battery that are damaged.

Not make a short circuit in the battery. A short-circuit can take place accidental when an object metallic (currency, clip or pen) cause one direct connection of the poles + and. of the battery (you throw metallic of the battery), for example when it takes a reserve battery in the pocket.

To make a short-circuit of the poles can damage the battery or the object that it is connected.

Not throw the batteries to the fire.

Throw the batteries following the local (for example, reciclaje) regulation.

Not throw the batteries like domestic garbage.

Avoid to load the battery under conditions of lack of air.

In order to maximize the yield of their battery:

Always use batteries and adapter of CA originals. The guarantee of the scale doesn't cover the damages caused as consequence of the use of other batteries and/or different loaders.

The tension of exit of adapting CA is of 9V, but the range of normal tension will oscillate between the 11V and the 15V. The new batteries or those batteries that have been stored during long periods of time can require a bigger time of load.

Maintain the battery to ambient temperature or a next temperature when it makes the load.

Not expose the batteries to inferior temperatures to -10°C or superiors at 45°C .

In the course of periods of long time, the batteries go losing load capacity progressively and they require longer times of recharge.

This is normal. If you load the battery regularly and he/she observes that it diminishes the period of operation or that it increases the period of load, it is probably moment of acquiring a new battery.

7. RS232 OUPUT

The BAR series scales can be ordered with an optional RS232 output.

Specifications:

- RS232 output of weighing data
- ASCII code
- 9600 Baud
- 8 data bits
- No Parity

Connector: 9 pin d-subminiature socket

- Pin 2: Input
- Pin 3: Output
- Pin 5: Signal Ground

Continuous output protocol.

		,			-/									k	g	CR	LF
-HEADER 1-			-HEADER 2-			-WEIGHT DATA-						-WEIGHT UNIT-		TERMIANTOR			

HEADER1: ST= STABLE, US= UNSTABLE

HEADER2: NT= NET, GS = GROSS

8. PARAMETERS

The scale has 9 parameters that can be set by the user plus a method of entering the calibration section. To set parameters press the  key during self checking.

The display will show the first function, "F1 UNIT". Pressing the  key will cycle through the other functions.

Pressing  key will allow you to set the function. Use the  key to leave a parameter unchanged.

Note: when display shows "TECH", press  key to enter, display shows "PIN", press , ,  key, and press  key to enter technical parameters setting mode.

PARAMETERS TABLE

FUNCTION	DESCRIPTION
F1 Unt	Sets the which units want to use: g / ct / lb / oz / d / gn / ozt / dwt / mm / tl.T / tl.c / tl.t / t use  to switch, use  key to sure.
F2 EL	Set the backlight to be on, automatic or off. EL on: Backlight always on EL Au: Backlight automatically turn on EL off: Backlight allways off Use  to switch, use  key to sure.
F3 Con	Set RS232 At first, select you will use RS232 or USB. \$ 232: use RS232 interface \$ USB: use USB interface Then will select interface work mode P Prt connect with mini printer P Cont connect with remote display or PC P Auto connect with mini printer, auto printer, auto print (auto accumulation) wirELE: wireless, when you use RW-W with kit I or BW-W with kit I, you can plug wireless kit I to balance RS232 interface, balance will send data to RW-W or BW-W interface. After you select communication mode, display will show b :xxx, this is baud rate, you can select from 600, 1200, 2400, 4800 or 9600 bps. If you select P Prt or P Auto (work with printer), after set baud rate, display will show tP or LP-50, this is printer mode, you can select normal mini printer (tP) or label printer (LP-50). If you select LP-50, display will show EnG or CH, this is print out language, you can select Chinese (CH) or English (EnG). Use  to switch, use  key to sure.
tech	Pin Enter the programming and calibration menus by entering the correct password ( ,  and ) and press  key to enter.
P1 Lin	Liner calibrate. After you press  key, display will show Pin again, press  ,  ,  key and  key to enter, display will show Load 0, remove all weight from the pan and wait until display show next step. Then display will show Load I, put the weight follow the display, wait for next step. Follow the display to load calibration weight on the pan. Until display show Load F, load full capacity weight on the pan. Note: calibration weight regulation: 0<load 1<load 2<load 3<load F(full capacity) Load1, Load2, Load3 need be integer weight.

P2 CAL	<p>Normal calibrate</p> <p>After enter cal mode, the display will show UnLoAd. Remove all weight from the pan and wait for stable.</p> <p>After stable, balance will show LoAd, put calibrate weight on the pan, here is the useable weight list.</p> <p>300g: 30/40/50/60/70/80/90/100/110/120/130/140/150/160/170/180/190/200/210/220/230/240/250/260/270/280/290/300</p> <p>600g: 60/80/100/120/140/160/180/200/220/240/260/280/300/320/340/360/380/400/420/440/460/480/500/520/540/560/580/600</p> <p>1500g: 150/200/250/300/350/400/450/500/550/600/650/700/750/800/850/900/950/1000/1050/1100/1150/1200/1250/1300/1350/1400/1450/1500</p> <p>3000g: 300/400/500/600/700/800/900/1000/1100/1200/1300/1400/1500/1600/1700/1800/1900/2000/2100/2200/2300/2400/2500/2600/2700/2800/2900/3000</p> <p>after stable, dosplay will show pass, calbrate complete.</p> <p>No need press any jey during calibrate.</p>
P3 Cnt	This display will show xxxxxxx for indicating the internal counts, you can press 0 key to escape.
P4 A2n	This option is used to select the auto zero maintain. Options: OFF, 0.5d, 1d, 2d, 4d, Use MODE key to shift, press 0 key to sure, press 0 key to escape.
P5 GrA	<p>This is used to adjust gravity, when move to another place with big latitude change, you can use this function to adjust.</p> <p>After press 0 key to enter this function, you can use 0, %, 0 key to key in the latitude rate (use 0 and % key move active digits, use 0 key increase value), after input, press 0 key to sure.</p> <p>Latitude rate= local latitude/the latitude of the calibrate place.</p> <p>For example, local latitude is 9.85, the latitude of the calibrate place is 9.75, by this way, the reading data will be bigger than standard, you need input 1.01025, then reading will be correct.</p> <p>Note: after you do calibrate, this rate will reset to 1.00000.</p>
P6 CAP	Set capacity, press MODE key to select: 300g/600g/1500g/3000g/6000g, for high resolution model, you can select 150g/200g/1500g/2000g, after you select, press 0 to confirm.

9. APPENDIX

9.1. ERROR CODES

During the initial power-on testing or during operation it is possible the scale may show an error message. The meaning of the error messages is described below.

If an error message is shown repeat the procedure that caused the message, turning the balance on, calibration or other functions. If the error message still is shown then contact your dealer for further support.

ERROR CODE	DESCRIPTION	POSSIBLES CAUSES
Err 3	Calibrate weight error	Use incorrect calibrate weight
Err 4	Initial Zero is greater than allowed (4% of maximum capacity) when power is turned on or when the 0 key is pressed.	Weight on the pan when turning the scale on. Excessive weight on the pan when zeroing the scale. Improper calibration of the scale. Damaged load cell. Damaged electronics.
Err 6	A/D count is not correct when turning the scale on	Platform not installed. Load cell damaged. Electronics damaged.
Err 7	Percentage setting error	0.01% weight must > 0.5d

Err 8	external calibration weight error	More than 4% error as linear calibration
-------	-----------------------------------	--

9.2 WEIGHING UNIT LISTING

Nº	UNIT SYMBOL	NAME OF THE UNIT	CONVERSION
1	g	gram	
2	g	gram (high resolution)	
3	ct	carat net	= 0,2 g
4	lb	pound	= 453,59237 g
5	oz	ounce	= 28, 349523000 g
6	d	drachm	= 1,7718450 g
7	ozt	ounce troy	= 31,1034768 g
8	dwt	scruples	= 1,55517400 g
9	mom	momme	= 3,749996 g
10	tL.H	tael twn	= 37,49996 g
11	tl.c	tael chn	= 37,799375 g
12	tl.t	troy tael	= 37,4290018 g
13	t	tola	= 11,6638039 g
14	bt	bangladeshi tola	1bt= 11,6638039 g 1bt= 16Ana= 96Roti
15	n	newton	1n= 102,04 g

10. TROUBLE SELF CHECKING

PROBLEM	POSSIBLE CAUSE	COMMON SOLUTIONS
Display is blank, No turn on test	Mains power is turned off. Power supply not plugged in. Power supply faulty. Internal battery not charged. Display turned off.	Check power is going to the scale and switch is on. Verify the voltage going to the scale matches the power supply labels on the power module or scale.
Display blank after turn on test, Error message, Display is locked	Pan not installed. Unstable weight. Load cell damaged. Mechanics damaged.	Check the pans are installed correctly. Try turning the scale on again.
OL or ---- appears on display	Maximum capacity exceeded. Load cell or mechanics damaged. Power supply faulty.	Check the pans are installed correctly. Try turning the scale on again.
---- or NULL displayed	Weight on scale below permissible limit. Pan has been removed. Pan support not seated properly. Power supply faulty. Load cell or mechanics damaged.	Check the pans are installed correctly. Try pressing  key. Try turning the scale on again.
Display is unstable	Drafts or air currents. Obstruction under pan. Sample is moving (animal weighing). Vibrations through table. Temperature changed dramatically. Power supply faulty.	Verify the scale is in a acceptable location and on a good table. Verify the power supply is correct for the scale.

Weight value incorrect	Calibration error, Recalibrate. Linearity error, set Linearity. Unit calibrated with inaccurate weight. Balance not level. Obstruction between sample and cover. Wrong unit of weight displayed.	Calibrate again, paying special attention to the mass used, the stability of the scale, and the weighing units required. If linearity can be set by the user it will be described in the user manual. Check pan is installed correctly. Verify installation is acceptable.
Cannot use Full Capacity	Over load Stops hitting pan support or hitting bottom of load cell. Shipping screw not removed if applicable. Electronic problem on A/D. Parameters set incorrectly. Load cell or mechanics damaged.	Look for obstruction under pan, shipping screws, and verify pan installation. Check the weighing units used.
Not linear	Overload stops hitting too soon. Load cell or mechanics damaged. A/D damaged.	Look for obstruction under pan, shipping screws, and verify pan installation. If linearity can be set by the user it will be described in the user manual.
Off center loading error	Adjust mechanics. Overload stops not correct. Load cell damaged.	Look for obstruction under pan, shipping screws, and verify pan installation.
Battery will not charge	Mains voltage not present or too low. Charging circuit failure. Battery failure.	Verify the batteries are rechargeable types. Check power supply voltage is correct.

11. WARRANTY

This scale is warranted against defects of manufacturing and materials for a period of 1 year, from the delivery date.

During this period, GIROPÈS SL will take charge of repairing the scale.

This warranty does not cover defects or damaged caused by misuse, overloading or improper installation contrary to the recommendations described in this manual.

This warranty does not cover shipping costs for the reparation of the balance.

**DECLARACIÓN DE CONFORMIDAD
DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITE**

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BAXTRAN S.L.
Pol. Empordà International
C/F. Parcela 15-16
E-17469 VILAMALLA (Girona)
- SPAIN -

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To which declaration refers to, conform with the followings standards or other normative documents:

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89/336/CE directive on the electromagnetic compatibility.

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And the following harmonized norms/ normes harmonisées:

EN55022 Clase B

EN61000-4-2

EN61000-4-3

EN61000-4-4

EN60950

Jordi ribalta

Director General
General Manager

marca propiedad de /
marque propriété de /
trade mark propriety of:

Pol. Empordà Internacional
Calle F. Parcela 15-16
17469 VILAMALLA
(Girona) SPAIN
T. (34) 972 527 212
F. (34) 972 527 211